RNA isolation protocol

Protocol edited by Nikita Jhaveri, 8th May 2025

(Original protocols: Erik Andersen 2010, Gaotian Zhang 2023)

Wear gloves and use RNase-free reagents and plastic Wipe down your bench and pipettes with 70% ethanol before proceeding Use filter tips for RNA isolation

Collection of worms

- 1. Wash the animals off the plates with M9 in a clean, labeled 15 mL conical tube.
- 2. Centrifuge at 254 g for one minute (Andersen lab: Eppendorf 5810R, 1100 rpm).
- 3. Discard the supernatant and add ~ 10 mL M9. Repeat step two.
- 4. Discard the supernatant. Add 1 mL of M9, mix well by pipetting up and down 5 times, and transfer the animals to a clean, labeled microcentrifuge tube.
- 5. Centrifuge in the tabletop centrifuge (Andersen lab: Eppendorf 5417C) at 10,000 rpm for one minute.
- 6. Discard the supernatant.
- 7. Add four volumes (of the worm pellet) of Trizol (Invitrogen, catalog # 15596018).
- 8. Vortex vigorously until completely resuspended.
- 9. Store at -80°C until all samples are ready to be processed.

RNA isolation

- 10. Thaw the samples from -80° C at 37° C.
- 11. Flash freeze the slurry in liquid nitrogen or dry ice.
- 12. Thaw at 37°C and vortex for 30 seconds.
- 13. Repeat steps 11 and 12.
- 14. Add two volumes (of the starting worm pellet) of Trizol and vortex vigorously.
- 15. Add 100 µL of prepped sand*.
- 16. Vortex vigorously for 10 minutes at room temperature (22°C, 1500 rpm).
- 17. Add 0.2 mL of chloroform (Fischer, catalog # C298-500) and vortex for 3 minutes at room temperature (22°C, 1500 rpm).
- 18. Centrifuge (Andersen lab: Eppendorf 5417C) at full speed (14,000 rpm) for 3 minutes.
- 19. Transfer the top aqueous layer to a new labeled microcentrifuge tube, taking care not to disturb the interphase.
- 20. Add 0.5 mL of isopropanol and mix the tube well by short vortex.
- 21. Incubate the tube at room temperature for 8 minutes.
- 22. Transfer the tube to ice for 2 minutes.
- 23. Centrifuge (Andersen lab: Eppendorf 5417C) at full speed (14,000 rpm) for 10 minutes.
- 24. Remove the supernatant and add 1 mL of 75% ethanol (made with RNase-free water).
- 25. Vortex vigorously and spin at full speed (Andersen lab: Eppendorf 5417C) (14,000 rpm) for 3 minutes.
- 26. Remove supernatant.
- 27. Centrifuge again at full speed (Andersen lab: Eppendorf 5417C) for 30 seconds to spin down the residual wash. Pipette to remove the residual wash. Be careful not to disturb the RNA pellet.

- 28. Air dry for 3 minutes or until the pellet appears almost completely dry.
- 29. Resuspend in 50 µL RNase-free water. Make sure RNA is fully suspended.
- 30. Aliquot 15 μ L to a separate tube for QC and store on ice. Transfer the master RNA tube to the -80°C freezer.
- 31. Assay RNA concentration using 10 µL on the Qubit with Qubit[™] RNA XR Assay Kit (cat. Q33224, Invitrogen via Life Technologies).
- 32. Dilute the remaining 5 μl to a concentration of 50 ng/μL to 500 ng/μL, for quality and integrity of the assay by Bioanalyzer and Nano chip. Each chip can measure 11 user samples.

*Sand is from Sigma (#274739). To prep sand, wash 2x in 1 M HCl, wash ~8x in RNAse-free water (until pH is ~7.0), bake to dry in an 80°C oven for 2 hours or more. Pre-dispense ~100 μ L aliquots before starting RNA prep so that they are handy for quick addition to each sample.